

VetScan

telemedicine

TeleRadiology

Table of contents:

Page 1:

Contents and contact information

Page 2:

Information needed for a TeleRadiology request

Page 2 - 5:

User manual

Prices for services will be sent out to all connected clinics once a year, if you are missing the pricelist, please send us an email at mail@vetscan.com.

Please contact us at mail@vetscan.com
or call +47 48325000
if you have any questions

- Information needed for a TeleRadiology request:

- Owner information:
 - First and last name
- Patient information:
 - Name
 - Sex
 - Species
 - Breed
 - Age
- Write all relevant information in “History”
- Follow the instructions in the user manual.

The screenshot shows a web interface for 'Nexus | Cases' with a 'TeleRadiology' section. It includes a 'Created on' date of 'Nov 18, 2019 - 11:56 AM', a 'Routine' dropdown menu, and a 'Creation in progress' status. Below this is a 'Patient' section with a search icon and input fields for 'Last Name', 'First name', 'Date of birth' (with a 'mm/dd/yyyy' placeholder), 'Sex' (with a dropdown arrow), and 'ID'.

- User manual:

The image displays two documents side-by-side. The left document is a blue 'User manual for TeleRadiology' with white text. It instructs users to 'Send X-Ray, CT or MRI images for evaluation' and provides 'Step 1: Log in at telemed.vetscan.com with the clinics username and password.' It also includes a '#' symbol and contact information: 'If you are missing username/password please email us at mail@vetscan.com'. The right document is a white login page for 'etiam One-Click Telemedicine' and 'AniCura telemedicine'. It features a 'Telemedicine Services' header, a link for 'New user AniCura Telemedicine, click here', and another link for 'Click here for user manual.' Below are 'Login' and 'Password' input fields and a 'Log in' button. The background of the right document shows a blurred image of a woman in a white lab coat.

STEP 2: 3 ways to add DICOM images, you can choose one of them (A, B or C):

Alternative A: Send directly from your PACS/X-Ray system.
If your clinic is not connected, please contact us and we will help you get started.

When the images are sent you will find them under «Exams»


Choose one or more exams and then press «create a case»
Choose TeleRadiology

Tip:
If you have more than one series of images, or old images that you want to compare, you can add them in the same case.


Alternative B: DRAG/DROP or from browser.
Save your DICOM images on your computer or other storage.
Go to «exams». Choose «Add exams»

Make sure you are in the area «From a directory»
Drag your DICOM images into the grey area or find them in a browser.
Press create a case – choose TeleRadiology.
You can only add DICOM images. Other files can be added at «Attached documents».
Please do not add X-Rays in other format than DICOM.


Alternative C: Add DICOM in Cases.
 Make sure to be in the page «cases».
 Press create a case – choose TeleRadiology



Go to the area DICOM exams and press «Add»
 Make sure you are in the area «From a Directory»
 Drag/drop you Dicam images to the grey area. Press «Add»



Tip: You can only add DICOM images. Other files can be added at
 «Attached documents».
 Please do not add X-Rays in other format than DICOM.



Tip:
 You will find a link to information
 and user manual above «History»
 #
 Remember to fill inn/check breed, age,
 ID and name.
 #
 If You dont have time to finish your case,
 the draft will save automatically.

- # You will get the report within 24 Hours (STAT 6 Hours).
- # The report will be sent to the e-mail registered in preferences.
- # You can change or add email in top right corner at «My Preferences»
- # You will also find the report and articles at telemed.vetscan.com
- # You can add follow-up questions when the report is finished, follow the instructions below.

etiam Nexus | Cases

Creation date Case author Case type Find test 1

To	State	Prio	Patient's name	Date of birth	Case type	Date updated	Case author
<input checked="" type="checkbox"/>	Report received, Report s...		Test 1 cardiology	Nov 1, 2018	Teleradiology	Nov 28, 2018	Clinical Case Manager Silje Ballestad Moland (...)
<input type="checkbox"/>	Taken over		test 2 cardiology	Nov 1, 2018	Teleradiology	Nov 28, 2018	Clinical Case Manager Silje Ballestad Moland (...)

How to add Questions:
 1: Make sure you are in «cases» in top left corner.
 2: Find your patient and mark it. Press open.

TeleRadiology

Created on Nov 28, 2018 - 2:00 PM

Routine Report received

Test 1

Date of birth: Nov 1, 2018
 Sex: M
 ID: 123

❖ You can add questions in all finished reports.
 ❖ You will get the answer by email, and it will be saved in the case.
 ❖ You can add documents together with the question if needed.

Contacts: Clinical Case Manager Silje Ballestad Moland

+ Add a step

4: Press «add a step» and then «follow-up question».

Follow-up question

Request taken over

Report

Follow-up question

5: Write your question in english.

Question: Write your Question.

Attached documents

6: Press send.

Send Remove

VetScan

telemedicine