

VetScan

telemedicine

TeleNeurology

Table of contents:

Page 1:

Contents and contact information

Page 2:

Information needed for a TeleNeurology request

Page 3 - 6:

User manual

Prices for services will be sent out to all connected clinics once a year, if you are missing the pricelist, please send us an email at mail@vetscan.com.

Please contact us at mail@vetscan.com
or call +47 48325000
if you have any questions

- Information needed for a TeleNeurology request:

- Owner information requested:
 - First and last name
- Patient information requested:
 - Name
 - Sex
 - Species
 - Breed
 - Age
- Choose your request.
 - Full clinical case. Turnaround time 24 hours.
- Write all relevant information in “History”
- Follow the instructions in the user manual.

The screenshot shows a web interface for 'Nexus | Cases'. At the top, it says 'TeleNeurology' and 'Created on Nov 18, 2019 - 11:56 AM'. Below this, there are two dropdown menus: 'Routine' and 'Creation in progress'. A section titled 'Patient' contains several input fields: 'Last Name', 'First name', 'Date of birth' (with a placeholder 'mm/dd/yyyy'), 'Sex' (a dropdown menu), and 'ID'.

- User manual:

The image shows two overlapping panels. The left panel is a blue box with white text titled 'User Manual for TeleNeurology'. It contains the following text: 'Send X-Ray, CT, MRI and videos for evaluation.', 'Step 1: Log in at telemet.vetscan.com with the clinics username and password.', a hash symbol '#', and 'If you are missing username/password please email us at mail@vetscan.com'. The right panel is a white login form for 'etiam One-Click Telemedicine'. It features the 'etiam' logo, the text 'Telemedicine Services', a link for new users, and a link for the user manual. Below this are input fields for 'Login' and 'Password', and a 'Log in' button. The background of the entire image shows a blurred photo of a woman and a man in a clinical setting.

STEP 2: 3 ways to send DICOM images, you can choose one of them (A, B or C):

Alternative A: Send directly from your PACS/X-Ray system.

If your clinic is not connected, please contact us and we will help you get started. When the images are sent you will find them under «Exams»

Choose one or more exams and then press «create a case»
Choose TeleNeurology

Tip:
If you have more than one series of images, or old images that you want to compare, you can add them in the same case.

Out:	Patient's name	Date of birth	Exam date	Modality	Description
<input checked="" type="checkbox"/>	test	Jul 15, 2015	Apr 23, 9:15 ...	CR	
<input type="checkbox"/>	test		Feb 18, 10:37...	CR	
<input type="checkbox"/>			Feb 1, 2013	CR	Tibia fibula,Abdomen,Thorax
<input type="checkbox"/>	test		Nov 21, 2018	DX	

Alternative B: DRAG/DROP
Save your DICOM images on your PC or other storage.
Go to «exams». Choose «Add exams»

Click an exam to show its preview...

Out:	Patient's name	Date of birth	Exam date	Modality	Description
<input type="checkbox"/>	test	Jul 15, 2015	Apr 23, 9:15 ...	CR	
<input type="checkbox"/>	test		Feb 18, 10:37...	CR	

Make sure you are in the area «From a directory»
Drag your DICOM images into the grey area. It may take a few minutes.
Press create a case – choose TeleNeurology.


Drop the files here

Patient	Date of birth	Sex	Patient ID	Exam date	Modality	Exam description	No. of files	Size


Community VetScan

Add Create a case Send to workstation/PACS Cancel


Alternative C: Add DICOM in Cases.
Make sure to be in the page «cases».
Press create a case – choose TeleNeurology



Go to the area DICOM exams and press «Add» Make sure you are in the area «From a Directory»
Drag your DICOM images into the grey area. It may take a few minuts.
Press "Add"



Tip: You can only add DICOM images. Other files can be added at «Attached documents».
Please do not add X-Rays in other format than DICOM.



Turnaround time: 24 Hours.

Remember to fill in/check breed, age, ID
and name.

If you dont have time to finish your case,
the draft saves automatically.

- # You will get the report within 24 Hours.
- # The report will be sent to the e-mail registered in preferences.
- # You can change or add email in top right corner at «My Preferences»
- # You will also find the report and eventual articles at telemed.vetscan.com
- # You can add follow-up questions when the report is finished, follow the instructions below.

etiam Nexus | Cases

272 Silje Ballestad Moland Clinical Case Manager

Creation date Case author Case type Find test 1

To	State	Prio	Patient's name	Date of birth	Case type	Date updated	Case author
<input checked="" type="checkbox"/>	Report received, Report s...		Test	Nov 1, 2018		Nov 28, 2018	Clinical Case Manager Silje Ballestad Moland (...)
<input type="checkbox"/>	Taken over		test 2	Nov 1, 2018		Nov 28, 2018	Clinical Case Manager Silje Ballestad Moland (...)

How to add Questions:

- 1: Make sure you are in «cases» in top left corner.
- 2: Find your patient and mark it. Press open.

TeleNeurology

Created on Nov 28, 2018 - 2:00 PM

Routine Report received

Test 1

Date of birth: Nov 1, 2018
Sex: M
ID: 123

Contacts

Clinical Case Manager Silje Ballestad Moland
VetScan - VetScan Telemed

4: Press «add a step» and then «follow-up question.»

Follow-up question Clinical Case Manager Silje Ballestad Moland (Clinica... Nov 28, 2018 - 2:00 PM

Request taken over Nov 28, 2018 - 9:02 PM

Report Nov 28, 2018 - 9:20 PM

Follow-up question Draft saved automatically

Question: Write your Question.

5: Write your question in english.

Attached documents

Add...

Draft saved automatically

6: Press send. Send Remove

- ❖ You can add questions in all finished reports.
- ❖ You will get the answer by email, and it will be saved in the case.
- ❖ You can add documents together with the question if needed.

VetScan

telemedicine